

TAKTVOLL

Intelligent Manufacturing in China



Electrosurgical Generator For Veterinary Use

ES-100V

Small Volume And Animal Friendly
Small size, easy to carry, cost-effective

Electrosurgical Generator

ES-100V For Veterinary Use

Integrated 6 different monopolar and bipolar modes and packed with dependable safety features, the ES-100-V satisfies the veterinarian's demands with **precision, safety, and reliability**.



ES-100V Feature

3 Cut Mode

PURE

Provides a standard Cutting effect.

BLEND 1

Provides a combination of 70% Cutting and 30% Coagulation

BLEND 2

Provides a combination of 50% Cutting and 50% Coagulation

2 Coag Mode

SPRAY

Designed to provide hemostasis in a broader area with less depth of tissue damage when compared to desiccation.

FORCED

Fast and effective standard coagulation with moderate thermal penetration. Slight carbonization effects.

1 Bipolar Mode

Standard

Bipolar coagulation which gently coagulates tissue for minimal tissue damage reducing carbonization and sticking.

TAKTVOLL[®]

Beijing Taktvoll Technology Co., Ltd.

Address: Room 10a-302, No. 13, South Fourth Street, Jingsheng, Jinqiao Science And Technology Industry Base, Tongzhou District, Beijing, China.

Tel: +86(10) 56495579 / 56495580 Web: www.taktvoll.com.cn E-mail: Support@taktvoll.com.cn