

TAKTVOLL®

**NEW PRODUCT!** ES-100V Plus LCD  
Touchscreen Electrosurgical Generator  
**ES-100V Plus**




The ES-100V Plus satisfies the veterinarian's demands with precision, safety, and reliability.



# Electrosurgical Generator

## ES-100V Plus For Veterinary Use



- Glass full version design, which is extremely tense and beautiful
- Color touch screen operation panel, flexible and easy to operate
- 6 Working modes
- Large blood vessel sealing function
- Electrosurgical Pens and Foot Switch Control
- Recently Mode, Power and Other Parameters
- Quick Adjustment of Power
- Cut and Coag in an Intermittent Manner



## ES-100V Plus Feature

3 Cut Mode		2 Coag Mode	1 Bipolar
<b>PURE</b>	Provides a standard Cutting effect.	<b>SPRAY</b>	<b>Standard</b> Bipolar coagulation which gently coagulates tissue for minimal tissue damage reducing carbonization and sticking.
<b>BLEND 1</b>	Provides a combination of 70% Cutting and 30% Coagulation	<b>FORCED</b>	
<b>BLEND 2</b>	Provides a combination of 50% Cutting and 50% Coagulation		