

TAKTVOLL

Intelligent Manufacturing in China



ES-120LEEP

Electrosurgical Generator Advanced Unit In Gynecology

The ES-120LEEP Integrated System combines three components — ES-120LEEP Generator, ES-120LEEP Smoke Evacuator and the ES-120LEEP Cart.

ES-120LEEP Integrated System

- Smoke evacuator eliminates odor at 99.99 percent efficiency level
- Fully integrated LEEP station with storage
- Heavy-duty casters facilitate mobility
- Elegantly designed for the modern procedure room



ES-120LEEP Feature

4 Cut Mode

PURE

Provides a standard Cutting effect.

BLEND 1

Provides a combination of 70% Cutting and 30% Coagulation

BLEND 2

Provides a combination of 50% Cutting and 50% Coagulation

BLEND 3

Provides a combination of 30% Cutting and 70% Coagulation

2 Coag Mode

FORCED

Fast and effective standard coagulation with moderate thermal penetration. Slight carbonization effects.

SOFT

Gentle coagulation with deep penetration, without carbonization, resulting in minimum adhesion of the electrode. Supported by the power control.

2 Bipolar Mode

Standard

Bipolar coagulation which gently coagulates tissue for minimal tissue damage reducing carbonization and sticking.

Precise

For microsurgical coagulation in the lower power range. Precise power settings and effects.

TAKTVOLL[®]

Beijing Taktvoll Technology Co., Ltd.

Address: Room 10a-302, No. 13, South Fourth Street, Jingsheng, Jinqiao Science And Technology Industry Base, Tongzhou District, Beijing, China.

Tel: +86(10) 56495579 / 56495580 Web: www.taktvoll.com.cn E-mail: Support@taktvoll.com.cn