

# TAKTVOLL

Intelligent Manufacturing in China



## Vessel Sealing Generator For Veterinary Use

ES-200VL

Both affordable and portable, easy-to-use. The platform is intended for open and laparoscopic surgical procedures.

# Vessel Sealing Generator

## ES-200VL For Veterinary Use

- Multi-mode with Ligasure.
- Automatic shutdown when sealing is done, without operator' s involvement.
- Intelligent power regulation, with lower energy to achieve best sealing on target tissue.



## ES-200VL Feature

### 3 Cut Mode

#### PURE

Provides a standard Cutting effect.

#### BLEND 1

Provides a combination of 70% Cutting and 30% Coagulation

#### BLEND 2

Provides a combination of 50% Cutting and 50% Coagulation

### 3 Coag Mode

#### SPRAY

Designed to provide hemostasis in a broader area with less depth of tissue damage when compared to desiccation.

#### FORCED

Fast and effective standard coagulation with moderate thermal penetration. Slight carbonization effects.

#### SOFT

Gentle coagulation with deep penetration, without carbonization, resulting in minimum adhesion of the electrode. Supported by the power control.

### 2 Bipolar Mode

#### Vessel Sealing

Seal and cut vessels up to and including 7mm in diameter.

Sealing vessels and tissue bundles without dissection or isolation , easy to operate.

#### Precise

For microsurgical coagulation in the lower power range. Precise power settings and effects.

**TAKTVOLL**<sup>®</sup>  
ES-200VL

Beijing Taktvoll Technology Co., Ltd.

Address: Room 10a-302, No. 13, South Fourth Street, Jingsheng, Jinqiao Science And Technology Industry Base, Tongzhou District, Beijing, China.

Tel: +86(10) 56495579 / 56495580 Web: [www.taktvoll.com.cn](http://www.taktvoll.com.cn) E-mail: [Support@taktvoll.com.cn](mailto:Support@taktvoll.com.cn)