

TAKTVOLL®

NEW GENERATION CO2 FRACTIONAL LASER MACHINE

LASER-3000

4 IN 1 FULLY-INTEGRATED



www.taktvollmed.com

LASER-3000

The CO2 Fractional Laser is an advanced fractional CO2 skin peeling laser system with a 10600nm wavelength. Apart from its gentle skin-peeling effect, it can deeply penetrate the laser beam into the dermis. This system is highly beneficial for skin recovery and can lead to long-term collagen remodeling, as well as improvement in skin conditions caused by exposure to light.

German import scanning head

- Long lifespan
- Strong energy
- Deep penetration
- Accurate scanning
- Low failure rate

Trachea

- Improve treatment experience and safety
- which has the functions of cooling and dust removal

Super powerful laser rod

TAKTVOLL technology enables the deepest penetration with the lowest energy. This reduces patient discomfort, downtime and thermal damage. Also enables the shortest pulses per energy level further enhancing patient comfort.

Light Guiding Arm

- 360° Smooth Operation
- 7 joints articulated light guiding arm
- Greatly reducing energy loss
- Convenient and flexible in operation
- Excellent use feeling

System Features

- Intelligent Operating System
- Multi-touch screen
- Intuitive touch-screen interface
- Displays treatment geometries
- Variety of treatment options





Fractional System

- Wrinkle removal
- Freckle spot removal
- Stretch marks removal
- All acne scars removal
- Skin resurfacing



Surgical System

- Medical cutting use
- Warts removal
- Naevus removal
- Acrochordon removal
- Polypus removal
- Tumors removal



Vaginal System

- Vaginal Tighten
- Vaginal Moistening
- Vaginal whitening
- Beautify vagina
- Urinary incontinence



Vulvae System

- Improve vulva redness
- Vulva skin rejuvenation

4 IN 1 FULLY-INTEGRATED


THE FUNCTION OF SYSTEM

SPOT SHAPE


Size adjustable suitable for different treatment area.

7 scan shapes: square, rectangle, round, triangle, oval, 6 –diamond shape and line

6 scan modes: sequence, random, maximum distance A&B, equidistance and fast shot




square



rectangle



round



triangle



diamond



oval



line


CO2 Fractional Laser Machine

Model	LASER-3000
Nominal Wavelength	10600nm
Laser power	60W
Work mode	Fractional, Surgical, Vaginal, Vulvae
Pulse width	pulse mode: 0.1ms-100ms fractional mode: 0.1ms-10ms
Beam delivery	seven joint Articulated Arm
Power	110V/220 V~ 50Hz

TAKTVOLL


Beijing Taktvoll Technology Co., Ltd



 Tel:+86-10-56495579

 Email:support@taktvoll.com.cn

 www.taktvollmed.com

 Headquarters:Room 302, 3rd Floor Floor, Building 10, No. 13A, Jingsheng South Fourth Street, Jinqiao Science and Technology Industry Base, Tongzhou District, Beijing, China.