

# TAKTVOLL

Intelligent Manufacturing in China


**99.99%**  
**Purification**  
**Efficiency**

## SMOKE-VAC 2000 Smoke Evacuator System



The Smoke-Vac 2000 medical smoking device adopts a 200W smoking motor to remove the harmful smoke generated effectively during gynecological LEEP, microwave treatment, CO2 laser and other operations.


## Smoke Evacuator System




## Accessories



Filter



Plastic Pipes



Sensor

## Specifications

Dimensions	26cm x 28cm x 12cm	Purification efficiency	99.99%
Weight	3.5kg	Degree of Particle Purification	0.3 $\mu$ m
Noise Levels	<60dB(A)	Operation Control	Manual/Auto/Foot Switch