

TAKTVOLL®

Intelligent Manufacturing in China


99.99%
Purification
Efficiency

SMOKE-VAC 2000 PLUS Smoke Evacuator System



The SMOKE-VAC 2000 PLUS features a new extended external filtration system, equipped with a 6.8-inch large elongated liquid crystal touchscreen, and introduces an added functionality of joint activation with ultrasonic scalpel systems.

Smoke Evacuator System



*Can connect two types of energy devices (high-frequency electrosurgical unit and ultrasonic scalpel) for simultaneous use during surgery. Filter core interface.

Accessories



Filter



medical smoke evacuation tubing



Sensor

Specifications

Dimensions	40cm x 39.5cm x 16cm	Purification efficiency	99.99%
Weight	8kg	Degree of Particle Purification	0.3µm
Noise Levels	<60dB(A)	Operation Control	Manual/Auto/Foot Switch/Electromagnetic