

# TAKTVOLL

Intelligent Manufacturing in China

99.9995%  
Purification  
Efficiency

## SMOKE-VAC 3000 Smoke Evacuator System



The new generation digital Smoke vac 3000 smoke evacuator system has low noise and strong suction. The turbocharging technology utilizes the suction power of the system, allowing you to perceive lightness, low noise, and effective smoke purification with a single touch.

# Smoke Evacuator System



## Accessories



## Specifications

Dimensions	355x197x248mm	Weight	7.3kg	Noise Levels	43.1 - 65.7dB
Flow Rate	1-3/8" (35MM) -76CFM	Degree of Particle Purification		0.12µm	
	1-1/4" (32MM) -74CFM	Operation Control		Manual/Auto/Foot Switch	
	7/8" (22MM) -38CFM	Suction Control		10%—100%	
	1/4" (6MM) -4.9CFM	Delay Time		0-60s	