


Model	SMOKE-VAC 2000 PLUS		
Description	Smoke evacuator (NEW PRODUCT!)		
Power Supply	100-240 VAC, 50/60 Hz		
Indications for Use	Suitable for use during open and laparoscopic surgeries to suction and purify smoke produced by energy devices including high-frequency electrosurgical units, laser equipment, and ultrasonic scalpels.		
Features	6.8-inch large elongated liquid crystal touchscreen		
	introduces an added functionality of joint activation with ultrasonic scalpel systems		
	Enhance High efficiency and durable filter		
	Filter monitoring function		
	Four working modes: handcontrol, foot control, auto, Sensor activation		
	Maximum output power: 500W		
	Rated airflow: 0–500 L/min		
	Filtration efficiency: 99.99% @ 0.3 μm		
	Noise level: 40–63 dB		
	Multiple adjustment of flow rate		
Standard	Standard configurations include the parts listed below.		
Model	Description	Quantity	Remark
SMOKE-VAC 2000 PLUS	Smoke Evacuation Unit	1pc	
ESU-510	Enhance Smoke filter	2pcs	
SJR-2553	filter TUBE, 200cm	1pc	

SJR-XYDB-001	electrosurgical pencil with smoke tube	1pc	
ESU-511	Laparoscopic Smoke Filtration tube	1pc	
SJR-256	filter adapter	1pc	
SJR-2734	Foot Switch	1pc	
SJR-644	Automatic Activation Device	1pc	
	AC Power cord	1pc	
	Instruction Manual	1pc	
	Certificate of quality	1pc	